



The everyday electronic magnifier

When your eyesight begins to fail and your glasses are no longer enough, reading can become a tremendous challenge. This is where the electronic magnifiers VISIO 24 come in, helping you retain that all important independence.

The VISIO 24 provide extra magnification with a crisp, clear picture so you can lead a normal life again. VISIO 24 magnify printed texts but also photos and other objects in full colour or contrasting colours if desired.

- Large 24" TFT widescreen monitor
- Crisp, clear picture at the press of a button
- 2MP / 1080p Full HD
- height-adjustable monitor
- Intuitively operated buttons with clear markings
- Reading line and reading curtain
- Easy to use – Switch-on and go
- Complete standard functionality
- Energy-saving LED lighting

VISIO 24

Retain your daily independence
with the VISIO 24

Place your documents and objects on the XY-table and VISIO 24 magnifies them on a widescreen monitor. Size, contrast and brightness can be adjusted to suit your needs. All the controls are mounted at the bottom of the screen for easy access, and are large and logical. For difficult reading material such as newspapers, the VISIO 24 can enhance the text making it clearer. VISIO 24 is the standard model with 6 buttons and offers more contrasting colour combinations. VISIO 24 is the perfect solution for independent living.



The aesthetically pleasing and space-saving design fits all environments.

Technical specifications

VISIO 24

- Control panel with 6 buttons
- 14 artificial colours with autofocus
- Full-HD Camera

Magnification range:

- 2,2 to 60 times magnification

Image:

- Full-colour, artificial colour combinations, inverse mode

Monitor:

- Full-colour, 24" (60,96 cm)
- 16:9 format, widescreen with LED backlight technology

Settings:

- Magnification, automatic contrast, image optimizing function
- Reading line and reading curtain

Voltage:

- 100 - 240 VAC 50 - 60 Hz 1.7 A

Size and weight:

- 53 cm x 63 cm x 58 cm
- appr. 11,5 kg (appr. 25.4 lbs.)



Easy to use function keys on the front.



Switch-on and VISIO 22 is immediately ready to use.

VisioBraille GmbH

Oßmaritzer Straße 4c phone.: 03641 - 2816-400
D-07745 Jena fax.: 03641 - 2816-416

e-mail: info@visiobraille.de
web: www.visiobraille.de

